iVision PROTECT

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iVision PROTECT

EYE DROPS with cross-linked
hyaluronic acid 0,2%
and PHMB 5 ppm

PROPERTIES: the product contains Hyaluronic Acid Cross-Linked and PHMB. Hyaluronic acid is a highly tolerated natural polymer, which is present both in lacrimal film and in human eye and has the property of soothing irritation of eye surface, while maintaining it optimally moistened. The Cross-linked formula improves the hyaluronic acid properties by achieving better bioavailability and resistance to degradation. iVision PROTECT also contains Polyhexamethylene Biguanide (PHMB), a substance that helps to inhibit microbial proliferation. iVision PROTECT is a lubricating, protective and moisturizing ophthalmic solution. Its use gives a ready feeling of relief, freshness and well-being.

USE INDICATIONS: iVision PROTECT eye drops are indicated to relieve eye redness, irritation, itch and dryness due to: air conditioning and heating, stay in smoky places or with other kinds of air pollution, prolonged exposure to wind or sun rays, chlorine in case of swimming practice in swimming pools, or saltiness if such activity is practiced in the sea, sport activities in general, deep concentration in study or work activities, use of contact lenses, post operating traumas, external bodies, excessive use of screen terminals and tiredness. iVision PROTECT is also indicated in case of blepharitis, trauma and in the post-operative course of ocular surface surgery (cataract, refractive, IVT, etc.).

USE INSTRUCTIONS: Put one or two drops in each eye squeezing the bottle gently. The product can be used every day, more times a day, according to necessities, also while wearing contact lenses. iVision PROTECT use

is simple, safe and without side effects; it can be used by everyone, including children and elderly.

PRECAUTIONS:

- The product is for external use only.
- The solution, once opened, must be used within 90 days.
- Do not use in case of damaged seal.
- Never put the dropper in contact with any kind of surface, including eye one, not to contaminate the product.
- Ensure the bottle is well closed after use.
- Store the bottle away from heat sources and in case under 30°C.
- Do not use after the expiry date printed on the packing.

WARNINGS:

- Do not use in case of specific hypersensitivity to some component. Keep far from children reach.
- In case of existing therapy with any ophthalmic drug seek professional advice before using iVision PROTECT contemporarily.
- Do not waste the bottle in the environment after use.

CONTRAINDICATIONS: There are no ascertained contraindications apart from allergy to specific components.

COMPLICATIONS: In case of trouble or persistent irritation, discontinue use and consult a specialist.

COMPOSITION: Polyhexamethylene Biguanide (PHMB) 0.0005%; Sodium Hyaluronate cross-link 0.2%; Hypromellose (HPMC) 0.2%; Disodium edetate, Borate buffer, Sodium Chloride, Excipients and purified Water up to 100%.

 Follow the instructions on the label strictly and keep this technical-informative note.

VALIDITY: 3 years if intact and properly preserved.