- TOALLITAS OFTÁLMICAS ESTÉRILES CHX
- LINGETTES OPHTALMIQUES STÉRILES CHX
- STERIELE OOGDOEKJES CHX
- STERILE OPHTHALMISCHE TÛCHER CHX
- TOALHITAS OFTÁLMICAS ESTÉREIS CHX
- SALVIETTE OFTALMISCHE CHX
- ΟΦΘΑΛΜΙΚΑ ΜΑΝΘΛΑΚΙΑ ΑΠΟΣΤΕΙΡΩΜΕΝΑ CHX

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Farmamix
Amposta, 20
Sant Cugat del Vallès
08174 (Barcelona) SPAIN

INDICATIONS
IVISION STERILE OPHTHALMIC WIPES CHX is indicated for topical treatment of the eyelids and eyelashes for both adults and children, even in the presence of secretions resulting from allergies, inflammation and/or infections, such as blepharitis and blepharoconjunctivitis of various kinds. Its use is particularly indicated for the pre-and/or post-ophthalmic surgery treatment. IVISION STERILE OPHTHALMIC WIPES CHX is also recommended for daily treatment of the periorcular area in individuals who have a predisposition for the occurrence of blepharitis and blepharoconjunctivitis of various kinds. Due to its peculiar composition, is also suitable for contact lens wearers.

PROPERTIES

All components of IVISION STERILE OPHTHALMIC WIPES CHX are useful to obtain optimal treatment, while fully respecting the particular physiology of the periorcular area, which is particularly delicate and sensitive. The product is capable of providing accurate treatment, by mechanical removal of secretions, foreign bodies, dust, etc.

COMPOSITION

Perforated wipes made of 100% cotton, wrapped in single dose twin-sachets and impregnated with 2,8 g of an aqueous solution: Aqua, Polysorbate 20, Sodium hyaluronate, Disodium EDTA, Trisodium citrate dihydrate, Aloe vera gel powder, Glycerol, Chlorhexidine digluconate, Citric acid mono-hydrate, Allantoin, Sodium hydroxymethylglycinate.

INSTRUCTIONS FOR USE

Having cleaned your hands carefully, open the sachet and extract the wipe. Apply the wipe to the eyelid area, wiping it delicately to allow for adequate cleansing, removal of ocular secretions, and any other types of deposit (smog, airborne particles, dust, etc.). In the case of particularly persistent secretions or the occlusion of a tear duct, it is advisable to warm the closed sachet for a few minutes in lukewarm water (approx. 40 °C) and then use the wipe as a hot compress, massaging the area delicately. On the other hand, in case of phlogosis, it is advisable to cool the sachet in the refrigerator for a few hours and then apply the open wipe to the eye, keeping it in place for a few minutes as a cold compress.

WARNINGS

Product sterility is guaranteed only if the sachet is intact. The wipe is disposable. Reusing it may spread infection. Keep in a cool and dry place. Keep out of reach of children. For external use. Discontinue

treatment in the case of evident irritation or hypersensitivity to the product. In case of infection do not use the same wipe for both eyes. Read instructions carefully.

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