DESCRIZIONE

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Rev. 01 del 30/11/2022

DESCRIPTION

repa**LysYal®** is a medical device coadjuvant in the healing of chronic wounds, ulcers and decubitus, due to the activation of granulation tissue that can promote the trophic process central to skin regeneration.

repa**LysYal®** is a cream for topical use that contains **T-LysYal®**, a supramolecular aggregate composed of Hyaluronic Acid, Lysine, Thymine and Sodium Chloride capable, with its hydrodynamic, of acting on the structures of the dermis and epithelium by modulating the inflammation process and the repair of dermo-epidermal lesions, preventing and obtaining a remodeled and less fibrotic scar tissue. The result of its great hydration is the ability to reboot the regenerative processes at the cellular level.

T-LysYal®, present in the device, increases tissue hydrodynamics, enhancing the hygroscopic capacity of the Hyaluronic Acid polymer chain, a fundamental structural component of the extra-cellular matrix (ECM). This helps to keep the matrix hydrated and supports its structural integrity as well. Given the high-level hygroscopicity, **T-LysYal®** solutions increase molecular transport capacity, closely reproducing the physiological action of endogenous hyaluronic acid, inducing increase in the hydration of the derma, facilitation of cell turnover and removal of free radicals; these properties facilitate the repair and regeneration processes of dermo-epidermal lesions.

INDICATIONS

Adjuvant in the repair process of skin lesions: accelerates the healing and regeneration of the skin, hydrates the superficial layers of the skin by restoring the cell-osmotic balance, preventing and remodeling the scar fibrotic tissue.

repa**LysYal®** facilitates the healing and regeneration process in cases of:

- Chronic skin lesions
- Radiotherapy lesions
- Pressure-sore lesions
- Diabetic ulcers
- Vascular ulcers
- Skin ulcers from autoimmune diseases (sclerodermas)
- Post-operative or post-trauma complications
- Burns

repa**LysYal®** is also effective in the presence of necrotic tissue through a hydrating and emollient action on the necrotic eschar. repaLysyal® facilitates the colliquation of the eschar thanks to its hydrolytic action, even in regions difficult to treat.

COMPOSITION

repa**LysYal®** cream: **T-LysYal®** 0,5% (Supramolecular aggregate of hyaluronic acid, lysine, thymine, sodium chloride) **Others ingredients:** aqua, ethylhexyl stearate, glyceryl stearate, dimethicone, propylene glycol, cetearyl alcohol, cyclopentasiloxane, glycolic acid, peg-100 stearate, triethanolamine, sorbitan stearate, stearic acid, methylparaben, propylparaben, sodium dehydroacetate, lecithin, tocopherol, ascorbyl palmitate, citric acid.

INSTRUCTION FOR USE

- Wash your hands before proceeding.
- Clean the wound with physiological solution.
- Apply a thin layer of the product on the entire affected surface, including the edges.
- Protect the area with a secondary dressing, preferably a greasy gauze where there is direct contact with the lesion.
- Repeat application once a day for the first week and then on alternate days.

WARNINGS – PRECAUTIONS FOR USE:

- The product is for external use only.
- Medical treatment with repa**LysYal®** calls for the respect of common hygiene and asepsis norms, therefore you must wash your hands carefully before proceeding with its application on injured skin.
- The affected area must only be washed using physiological solution in order to remove any foreign bodies, while excess secretion and/or blood must be stemmed with gauze compresses.
- Do not use in combination with other products
- Excess necrotic material and macerated edges of skin should be removed in turn by use of delicate curettage by a health professional.
- In presence of purulent material consult your doctor before proceeding with use.
- Be sure to cover the wound post application to prevent contact and external infections.
- Keep out of reach of children. Do not swallow.
- The expiry date indicated on the label refers to the closed and correctly stored product.
- Do not use after the expiry date. Do not use if the package is open or damaged.

CONTRAINDICATIONS AND SIDE EFFECTS

- This product is not suitable for subjects who are intolerant/hypersensitive to any of its ingredients.
- There are no known undesirable effects of clinical relevance.
- Immediately following application to injured skin, there may be a burning sensation that tends to fade by itself after a few minutes.
- Should there be undesirable effects other than those described, stop the application and consult your doctor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer SILDEHA and the competent authority of the Member State in which the user and/or patient is established.

PREGNANCY AND LACTATION

Although effects of repa**LysYal®** on pregnancy and lactation are not known, the use of the product in these conditions should be subordinated to the opinion of the physician.

STORAGE CONDITIONS


Keep at room temperature. Keep off heating sources. Microbiologically controlled device.

PACKAGING

Cream: tube 30 ml.
Cream: 6 tubes of 2 ml.

Rev. 01 del 30/11/2022

SILDEHA Pharma S.r.l. Via Egadi, 7 - 20144 Milano - IT
Products Licensed by SILDEHA Swiss SA (International patents)

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